

Requisition of HR-Urticaria Test



RefLab • P.O. box 590 • DK- 2200 Copenhagen N
Phone: +45 3545 7512 • Fax +45 3545 7506 • reflab@reflab.dk • www.reflab.dk

PATIENT	NAME OF HOSPITAL/CLINIC
<p>Name</p> <p>Address</p> <p>.....</p> <p>Date of birth</p> <p>Date of blood sampling</p> <p>Male <input type="checkbox"/> Female <input type="checkbox"/></p>	<p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>

Detection of autoantibodies against the IgE and/ or IgE receptor

HR-Urticaria Test \perp

Practical information – Please note

Serum (**1-2 mL per patient placed in a transport container**) is posted by ordinary mail at ambient temperature to RefLab.

Method: Patient serum and donor basophils are incubated. Histamine is measured in stimulated and un-stimulated cells and the results are expressed as histamine release in percentage of the cellular histamine content.

A histamine release > 16.5 % is a positive test result.

The test result will be returned to the clinic within 10 days after receipt to the laboratory.